

CATHTOP SOLUTION FOR DIAGNOSTIC





CATHTOP PUNCTURE NEEDLE



Features:

- 1. Sharp tip that has low resistance increases accessibility and reduce tissue trauma.
- 2. Short tip has better controllability.
- 3. Tapered hub inner diameter facilities insertion of the guide wire.
- 4. Needles of customized sizes are also available.

Cathtop Puncture Needle			
ORDERING INFORMATION			
LENGTH	DIAMETER	REFERENCE NUMBER	
7 CM	18G	CTPN-1807	
4 CM	20G	CTPN-2004	







2 - Port

Features:

- Available in 2 port Right ON/OFF and 3 port ON/OFF.
- 2. Withstands pressure up to 500 Psi.
- 3. Presence of rotation adapters.
- 4. Ease of inspection.



3-Port

Cathtop Manifold		
ORDERING INFORMATION		
PRODUCT VARIENT	REFERENCE NUMBER	
2 PORT RIGHT OFF	CTMF02RF	
2 PORT RIGHT ON	CTMF02RN	
3 PORT RIGHT OFF	CTMF03RF	

| CATHTOP | YCONNECTOR



Push-Pull Y connector

Features:

- 1. Two different Sideons of Y adapter suits various customers.
- 2. Push Pull Y connector makes it easier for the doctor to operate precisely.



Push-Click Y connector

Features:

- 1. Three-stage valve creates maneuverability without any bleeding.
- 2. Auto-closing mechanism ensures the device to move from semi-closed position to closed position automatically when contrast is injected from the sideon lumen.
- 3. Ergonomic design allows for one-hand operation.

Cathtop Y Connector		
ORDERING INFORMATION		
DESCRIPTION	REFERENCE NUMBER	
7F, Y CLICK, SIDEON TUBING AND STOPCOCK, INSERTION TOOL, RED/COPPER TORQUER, TYVEK-PLACTIC BLISTER	CTYC-PP-7F	
	CTYC-PC-7F	

CATHTOP INFLATION DEVICE



Normal Inflation
Device



Hybrid Inflation Device



Gun Type Inflation Device

Features:

- 1. Gun type inflation device which comes with a HP stopcock.
- 2. Comes with a 20ml clear syringe and 30atm/bars & a Precision manometer.
- 3. Ergonomic design to facilitate operation of the device.
- 4. Trigger release mechanism helps in rapid inflation and deflation.
- 5. Inclination of manometer provides easy and precise view of Pressure points.

Cathtop Inflation Device		
ORDERING INFORMATION		
DESCRIPTION	REFERENCE NUMBER	
25ML, 30ATM/BAR, HP-		
STOPCOCK, TYVEK-PLASTIC	CTID25-HT	
BLISTER		
20ML, 30ATM/BAR	CTID20-NT	
20ML, 30ATM/BAR WITH	CTID20-GT	
HP STOPCOCK		



